



Spirometrics



Now use Spirometrics D1010 Filters on other brand spirometers with the appropriate adaptor listed below!

Now is the perfect time to talk with your customers about their spirometry filter needs.

The Spirometrics PFT filter is considered one of the best on the market for reducing contamination and concern about the spread of infectious disease. It successfully blocks more than 99% of bacteria and viruses.*† Its new ergonomic design offers a better mouth fit for maximum test accuracy, meeting flow resistance requirements over the entire range of ATS flows from zero to 14 liters/sec.†



* Independently verified by Nelson Labs.
† For test data, please contact Spirometrics.

Filtration Efficiency
99.9% Bacterial*†
99.6% Viral*†

Cross Contamination Efficiency
99.9999% Bacterial
99.9984% Viral

Spirometrics Filter Adaptor

CAT#	FITS
A3001-2	SensorMedics 2100, 2200, 6200 Vmax, Pall size
A3002-2	Spirotech, Ohio, Collins Survey
A3003-2	Vitalograph
A3004-2	MultiSPIRO (older models), Morgan Benchmark DLC0, SensorMedics 2100, 2200, Schiller
A3005-2	Gould/SensorMedics (models 3400 & 2450), Welch Allyn
A3006-2	MultiSPIRO (current models), Creative Biomedics, Keystone (old), Cosmed pony, Vacumed pony
A3007-2	Cybermedic Spinaker Excel, MedGraphics Spiro Tube, Spirovit
A3008-2	Spirotech, MedGraphics 1085, Collins DS, GS
A3009-2	Cybermedic CM3, Gould, Jones Satellite
A3010-2	Riko, Accutrax
A3011-2	Tamarac, Burdick (old), CDX
A3013-2	Brentwood 4000
A3015-2	Adapts D1010 filter to calibration syringe
A3017-2	Collins GS, CPL
D1010 Filter is a direct replacement for the KoKo Moe™ filters, #810000 and #819000 and the V-MAX by Susquehanna Micro, #CYB-CB1503	
Fits all Koko spirometers (KoKo and KoKo-trek)	

Spirometrics

P.O. Box 680
Gray, ME 04039
1.800.767.0004
www.spirometrics.com